

15210

## Foam Anti-Fog Surgical Face Mask

- Complies with ASTM F2100 Level 1 Requirements
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Softer longer ties promote secure comfortable fit




### Specifications

Product	Foam Anti-Fog Surgical Face Mask
Color	Blue
Fluid Resistance	80 mm Hg
Filtration - BFE	≥98
Filtration - PFE	≥98
Delta P	≤2.8
Flame Spread	Class 1
Attachment Design	Vertical Ties
Material	Cellulose Inner and Polypropylene Outer
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Case

This product is NOT made with natural rubber latex

\*According to AAMI Standard PB70:2012 Liquid Barrier Performance



### Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 μm:	≥ 95%
PFE (%) @ 0.1 μm:	≥ 95%
ΔP (mm H <sub>2</sub> O/cm <sup>2</sup> ):	< 4
Flame Spread:	Class 1