

15930

## FluidGard® 120 Anti-Fog Surgical Mask with Extended Shield

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Softer longer ties promote secure comfortable fit



### Specifications

Product	Surgical Mask with Extended Shield
Color	Blue
Fluid Resistance	120 mm Hg
Filtration - BFE	≥99
Filtration - PFE	≥98.9
Delta P	≤3.3
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	25/Box, 4/Boxes, 100/Cs

This product is NOT made with natural rubber latex



### Complies with ASTM F2100

Fluid Resistance (mmHg):	120
BFE (%) @ 3.0 µm:	≥ 98%
PFE (%) @ 0.1 µm:	≥ 98%
ΔP (mm H <sub>2</sub> O/cm <sup>2</sup> ):	< 5
Flame Spread:	Class 1