



5030

Chemotherapy Gown

360° Full coverage.

- Tested for commonly used Chemotherapy drugs.*
- PE coated splash resistant front and sleeves.
- Polypropylene back panels help maintain coolness and comfort.
- Knit cuffs provide comfort.
- Secure hook and loop neck closure, and full back.



Specifications

Product	Chemotherapy Gown
Size	X-Large
Color	White
Wrist	Knit Cuffs
Material	PE Coated Polypropylene (front), Polypropylene (back)
Design	Hook & Loop Neck, Full Back
Seams	Sonically Sealed
Packaging	5/Bag, 6/Bags, 30/Case

This product is NOT made with natural rubber latex

*The Chemotherapy gown material has been tested against the following chemotherapy drugs and passes testing per ASTM D 6978-05 standards.

- Dacarbazine (DTIC)
- Cyclophosphamide (Cytoxan)
- 5-Fluorouracil
- Etoposide
- Paclitaxel (Taxol)
- Carmustine (BCNU)
- Doxorubicin
- Cisplatin (Platinol)
- Thio-TEPA

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.