

15130

KittyKare Procedure Face Mask

- Child friendly printed outer facing
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



Specifications

Product	KittyKare Procedure Face Mask
Color	KittyKare Print
Fluid Resistance	N/A
Filtration - BFE	≥99.1
Filtration - PFE	≥99.1
Delta P	≤3.1
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10/Boxes, 500/Cs

This product is NOT made with natural rubber latex