

15150

## Children's Face Mask

- Child-sized to improve fit for children
- Child friendly printed outer facing
- Soft comfortable earloops provide easier donning

### Specifications

Product	Children's Face Mask
Color	Smiley Face
Fluid Resistance	N/A
Filtration - BFE	≥99.1
Filtration - PFE	≥99.1
Delta P	≤3.1
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	75/Box, 10/Boxes, 750/Cs

This product is NOT made with natural rubber latex

