

15220

Lite Pouch Surgical Mask

- Pouch style mask provides off the face design giving more room to breathe
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit

Specifications

Product	Lite Pouch Surgical Mask
Color	Blue
Fluid Resistance	N/A
Filtration - BFE	≥99.8
Filtration - PFE	≥99.1
Delta P	≤2.9
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Cs

This product is NOT made with natural rubber latex

