

FSS Gown Products

Contract #

FSS: V797P-4398B

Expire

3/14/2021



Precept[®]

Optimizing Safety, Maximizing Comfort

Blue Film Gowns

Over-the-Head thumb loop gown with open back.

Item #	Description
51586	Blue Film Gown
51586B	<ul style="list-style-type: none">● Blue Polyethylene film* provides enhanced barrier protection.● Generous sizing improves comfort.● Over-the-Head donning improves compliance.● Tapered wrists and thumb loops ensure proper glove donning.● Thermal Impulse sealed seams.● Bulk packaging option available. 51586 - 15/Box, 75/Cs 51586B - Individually Bagged, 100/Cs● Size - X-Large

This product is NOT made with natural rubber latex

*Material Passes F1671 Viral Testing

ASTM F1671: "Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System"

For more information on this test method, please visit astm.org

