

15211

Film Anti-Fog Surgical Face Mask


- Complies with ASTM F2100 Level 1 Requirements
- Anti-fog film inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit



Specifications

Product	Film Anti-Fog Surgical Face Mask
Color	Green
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.7
Filtration - PFE	≥99.7
Delta P	≤2.9
Flame Spread	Class 1
Attachment Design	Soft Vertical Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Case

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 µm:	≥ 95%
PFE (%) @ 0.1 µm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.