

15320

## **FluidGard® 160 Anti-Fog Surgical Mask**


- Complies with ASTM F2100 Level 3 Requirements
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Softer longer ties promote secure comfortable fit



### Specifications

Product	FluidGard® 160 Anti-Fog Surgical Mask
Color	Blue Diamond
Fluid Resistance	160mm Hg
Filtration - BFE	≥99.8
Filtration - PFE	≥99
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Cs

This product is NOT made with natural rubber latex



### Complies with ASTM F2100

**Fluid Resistance (mmHg): 160**  
**BFE (%) @ 3.0 µm: ≥ 98%**  
**PFE (%) @ 0.1 µm: ≥ 98%**  
**ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5**  
**Flame Spread: Class 1**

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.