

14400

FluidGard® 160 Procedure Mask


- Complies with ASTM F2100 Level 3 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



Specifications

Product	FluidGard® 160 Procedure Mask
Color	Blue
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99
Filtration - PFE	≥99
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10 Boxes, 500/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 160
BFE (%) @ 3.0 µm: ≥ 98%
PFE (%) @ 0.1 µm: ≥ 98%
ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1