

15255

## FluidGard® 120 Horizontal Tie Surgical Mask with Anti-Glare Shield

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-fog, anti-glare treatment provides enhanced clarity
- Wider design provides more facial coverage than standard sized masks
- Soft longer horizontal ties promote secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



### Specifications

Product	FluidGard® 120 Horizontal Tie Surgical Mask with Anti-Glare Shield
Color	Blue with Pattern
Fluid Resistance	120 mm Hg
Filtration - BFE	≥99
Filtration - PFE	≥99
Delta P	≤3.3
Flammability	Class 1
Attachment Design	Horizontal Ties
Mask Construction	Sonic Welded
Packaging	25/Box, 4 Boxes, 100/Cs

This product is NOT made with natural rubber latex



### Complies with ASTM F2100

**Fluid Resistance (mmHg):** 120  
**BFE (%) @ 3.0 µm:** ≥ 98%  
**PFE (%) @ 0.1 µm:** ≥ 98%  
**ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>):** < 5  
**Flame Spread:** Class 1