

15311

FluidGard® 160 Anti-Fog Procedure Mask with Anti-Glare Shield


- Complies with ASTM F2100 Level 3 Requirements
- Coated film shield with proprietary anti-fog, anti-glare treatment provides enhanced clarity
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



Specifications

Product	FluidGard® 160 Anti-Fog Procedure Mask with Anti-Glare Shield
Color	Blue Diamond
Fluid Resistance	160mm Hg
Filtration - BFE	≥99.8
Filtration - PFE	≥99
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	25/Box, 4/Boxes, 100/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	160
BFE (%) @ 3.0 µm:	≥ 98%
PFE (%) @ 0.1 µm:	≥ 98%
ΔP (mm H ₂ O/cm ²):	< 5
Flame Spread:	Class 1