

15524

## FluidGard® 160 Laser Surgical Mask


- Complies with ASTM F2100 Level 3 Requirements
- > 99% efficiency filtering laser plume particles of 0.1 microns or larger
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit



### Specifications

Product	FluidGard® 160 Laser Surgical Mask
Color	Lavender
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99.7
Filtration - PFE	≥99.1
Delta P	≤3.3
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Cs

This product is NOT made with natural rubber latex



### Complies with ASTM F2100

Fluid Resistance (mmHg):	160
BFE (%) @ 3.0 µm:	≥ 98%
PFE (%) @ 0.1 µm:	≥ 98%
ΔP (mm H <sub>2</sub> O/cm <sup>2</sup> ):	< 5
Flame Spread:	Class 1