

15525

FluidGard® 160 Anti-Fog Laser Surgical Mask


- Complies with ASTM F2100 Level 3 Requirements
- > 99% efficiency filtering laser plume particles of 0.1 microns or larger
- Anti-fog foam inhibits eyewear fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit



Specifications

Product	FluidGard® 160 Anti-Fog Laser Surgical Mask
Color	Lavender
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99.7
Filtration - PFE	≥99.1
Delta P	≤3.3
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 160
BFE (%) @ 3.0 µm: ≥ 98%
PFE (%) @ 0.1 µm: ≥ 98%
ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1