

MSPV Products

Precept®

Optimizing Safety, Maximizing Comfort

Contract #

FSS: V797P-4398B

(MSPV) J & A Formulary

Expire

3/14/2021

Standard Surgical Face Masks

| Item # | Description |
|---------|---|
| 15220 | Lite Pouch Surgical <ul style="list-style-type: none">• Pouch style mask provides off the face design giving more room to breathe• Specially designed inner facing improves comfort• Integrated layers work together to maximize breathability• Soft longer ties promote secure comfortable fit• Blue, 50/Box, 6/Boxes, 300/Cs |
| 15250 | Lite & Cool Horizontal Tie Surgical Mask <ul style="list-style-type: none">• Wider design provides more facial coverage than standard sized masks• Soft longer horizontal ties promote secure comfortable fit• Specially designed inner facing improves comfort• Integrated layers work together to maximize breathability• Blue, 50/Box, 12/Boxes, 600/Cs |
| 65-3122 | Comfort-Plus™ Surgical Mask with Stretch Knit Ties <ul style="list-style-type: none">• Soft sewn stretch knit ties provide secure comfortable fit• Specially designed inner facing improves comfort• Integrated layers work together to maximize breathability• Green, 50/Box, 5/Boxes, 250/Cs |
| 65-3348 | DualGard™ Surgical Mask with Foam Shield <ul style="list-style-type: none">• Foam Shield™ Coated film shield with proprietary anti-fog treatment provides enhanced clarity• Soft sewn stretch knit ties provide secure comfortable fit• Specially designed inner facing improves comfort• Integrated layers work together to maximize breathability• Blue Diamond, 25/Box, 4/Boxes, 100/Cs |



15220



15250



65-3348

Standard Specialty Face Masks

| Item # | Description |
|---------|--|
| 65-3310 | Laser Mask with Stretch Knit Ties <ul style="list-style-type: none">• > 98% efficiency filtering laser plume particles of 0.1 microns or larger• Soft sewn stretch knit ties provide secure comfortable fit• Specially designed inner facing improves comfort• Integrated layers work together to maximize breathability• White, 50/Box, 5/Boxes, 250/Cs |
| 65-3322 | Tape Fog-Shield® Surgical Mask <ul style="list-style-type: none">• Fog-Shield® anti-fog tape inhibits eyewear fogging improving vision• Strong polyester ties provide secure comfortable fit• Specially designed inner facing improves comfort• Integrated layers work together to maximize breathability• Green, 50/Box, 5/Boxes, 250/Cs |



65-3322

MSPV Products

Precept[®]

Optimizing Safety, Maximizing Comfort

Contract #

FSS: V797P-4398B

(MSPV) J & A Formulary

Expire

3/14/2021

Face Shields

Item # Description

-
- | | |
|--------------|--|
| 15990 | 3/4 Length Face Shield <ul style="list-style-type: none">• Stretchable headband, with soft comfortable foam strip• Length 7 ½", 50/Cs |
| 15995 | Full Length Face Shield <ul style="list-style-type: none">• Stretchable headband, with soft comfortable foam strip• Length 9 ½", 50/Cs |
-



15990

Eyewear

Item # Description

-
- | | |
|--------------|--|
| 1599A | Eye Shield - Assembled <ul style="list-style-type: none">• Lightweight, wrap around eye shields provide protection from the front and sides• Eye shield lenses are optically clear and anti-fog• Can be worn over most prescription glasses• Convenient packaging• Assorted colors, 10/Bag, 5/Bags, 50/Cs |
|--------------|--|
-

