


15115

## Procedure Mask with Extended Shield


- Complies with ASTM F2100 Level 1 Requirements
- Coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



### Specifications

Product	Procedure Face Mask with Extended Shield
Color	Blue
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.6
Filtration - PFE	≥98.7
Delta P	≤2.4
Flame Spread	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	25/Box, 4/Boxes, 100/Case

This product is NOT made with natural rubber latex



### Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 µm:	≥ 95%
PFE (%) @ 0.1 µm:	≥ 95%
ΔP (mm H <sub>2</sub> O/cm <sup>2</sup> ):	< 4
Flame Spread:	Class 1