

15100

Isolation Face Mask

- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning

Specifications

Product	Isolation Face Mask
Color	Yellow
Fluid Resistance	N/A
Filtration - BFE	≥98.7
Filtration - PFE	≥98.9
Delta P	≤2.9
Flame Spread	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10/Boxes, 500/Cs

This product is NOT made with natural rubber latex

