

15121

Procedure Mask


- Complies with ASTM F2100 Level 1 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



Specifications

Product	Procedure Mask
Color	Yellow
Fluid Resistance	80 mm Hg
Filtration - BFE	≥98.4
Filtration - PFE	≥97.7
Delta P	≤3.5
Flammability	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10/Boxes, 500/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 μm:	≥ 95%
PFE (%) @ 0.1 μm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1