

15216

Sensitive Skin Foam Anti-Fog Surgical Mask


- Complies with ASTM F2100 Level 1 Requirements
- Color-free face mask with specially designed inner facing improves comfort
- Anti-fog foam inhibits eyewear fogging improving vision
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit



Specifications

Product	Sensitive Skin Foam Anti-Fog Surgical Mask
Color	White
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.7
Filtration - PFE	≥99.7
Delta P	≤2.9
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 6/Boxes, 300/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 µm:	≥ 95%
PFE (%) @ 0.1 µm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1