

15250

Lite & Cool Horizontal Tie Surgical Mask

- Wider design provides more facial coverage than standard sized masks
- Soft longer horizontal ties promote secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	Lite & Cool Horizontal Tie Surgical Mask
Color	Blue
Fluid Resistance	N/A
Filtration - BFE	≥99.7
Filtration - PFE	≥99.0
Delta P	≤2.1
Flammability	Class 1
Attachment Design	Horizontal Ties
Mask Construction	Sonic Welded
Packaging	50/Box, 12/Boxes, 600/Cs

This product is NOT made with natural rubber latex