

15331

FluidGard® 160 Anti-Fog Surgical Mask with Anti-Glare Shield


- Complies with ASTM F2100 Level 3 Requirements
- Coated film shield with proprietary anti-fog, anti-glare treatment provides enhanced clarity
- Anti-fog foam inhibits fogging improving vision
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Softer longer ties promote secure comfortable fit



Specifications

Product	FluidGard® 160 Anti-Fog Surgical Mask with Anti-Glare Shield
Color	Blue Diamond
Fluid Resistance	160mm Hg
Filtration - BFE	≥99.8
Filtration - PFE	≥99.0
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	25/Box, 4/Boxes, 100/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 160
BFE (%) @ 3.0 µm: ≥ 98%
PFE (%) @ 0.1 µm: ≥ 98%
ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1