

15830

FluidGard® 120 FluidPouch™ Surgical Mask with Anti-Glare Shield

- Complies with ASTM F2100 Level 2 Requirements
- Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity
- Pouch style mask provides off the face design giving more room to breathe
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit



Specifications

Product	FluidGard® 120 FluidPouch™ Surgical Mask with Anti-Glare Shield
Color	Blue with Pattern
Fluid Resistance	120 mm Hg
Filtration - BFE	≥99.1
Filtration - PFE	≥99.6
Delta P	≤2.4
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	25/Box, 4 Boxes, 100/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 120
BFE (%) @ 3.0 µm: ≥ 98%
PFE (%) @ 0.1 µm: ≥ 98%
ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1