

15900

FluidGard® 120 Procedure Mask

- Complies with ASTM F2100 Level 2 Requirements
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft comfortable earloops provide easier donning



Specifications

Product	FluidGard® 120 Procedure Mask
Color	Blue
Fluid Resistance	120 mm Hg
Filtration - BFE	≥99.0
Filtration - PFE	≥98.9
Delta P	≤3.3
Flame Spread	Class 1
Attachment Design	Earloop
Mask Construction	Sonic Welded
Packaging	50/Box, 10/Boxes, 500/Case

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	120
BFE (%) @ 3.0 µm:	≥ 98%
PFE (%) @ 0.1 µm:	≥ 98%
ΔP (mm H ₂ O/cm ²):	< 5
Flame Spread:	Class 1