

# Precept<sup>®</sup>

Optimizing Safety, Maximizing Comfort

MADE IN  
NORTH AMERICA

35-5703

## Tyvek<sup>®</sup> Gown

- Full back, generous sizing
- Sewn seams on sleeves, armholes, and sholders.
- Low-linting, anti-static.
- Secure neck ties and knit cuff.
- Lightweight Tyvek<sup>®</sup> material provides excellent protection against fluids and bacteria.
- Tyvek<sup>®</sup> is noted for its strength and resistance to tears and punctures.

### Specifications

Product	Tyvek <sup>®</sup> Gown
Size	Universal, one size fits most
Dimensions	Length 45", Width 56"
Color	White
Wrist	Knit Cuff
Material	Tyvek <sup>®</sup>
Design	Secure Neck Ties, Full Back
Seams	Sewn
Packaging	30/Case

This product is NOT made with natural rubber latex

Tyvek<sup>®</sup> is a registered mark of DuPont™

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.

