

# Precept®

Optimizing Safety, Maximizing Comfort

MADE IN  
NORTH AMERICA

AAMI  
Level  
2

51178

## Isolation Gown

Full Coverage  
AAMI Level 2

- Multi-layer SMS material\* delivers AAMI Level 2 protection and comfort.
- Secure neck ties.
- Elastic wrists.
- Sonic-welded seams.
- Full coverage with overlapping back panels.

### Specifications

Product	Isolation Gown, Full Coverage
Size	X-Large
Dimensions	Length 48", Width 60"
Color	Yellow
Wrist	Elastic
Material	Multi-Layer SMS
Design	Secure Neck Ties, Full Coverage
Seams	Sonic-Welded
Packaging	10/Bag, 10 Bags, 100/Case

This product is NOT made with natural rubber latex

**\*AATCC 42 - Water Resistance: Spray Impact Penetration Test**  
**AATCC 127 - Water Resistance: Hydrostatic Pressure Test**  
**According to AAMI Standard PB70:2012 Liquid Barrier Performance**

For more information on these test methods, please visit [aatcc.org](http://aatcc.org) and [aami.org](http://aami.org)

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.

