

Precept®

Optimizing Safety, Maximizing Comfort

AAMI
Level
2

MADE IN
NORTH AMERICA

59979

Isolation Gown

Full Coverage
AAMI Level 2 gown.

- Multi-layer SMS material* delivers AAMI Level 2 protection and comfort.
- Tape tab neck closure.
- Elastic wrist.
- Sonically sealed seams on sleeves, armholes, shoulders, and side ties.
- Full coverage with overlapping back panels.

Specifications

Product	Isolation Gown, Full Coverage
Size	Universal
Dimensions	Length 45", Width 55"
Color	Yellow
Wrist	Elastic
Material	Multi-layer SMS
Design	Tape Tab Neck Closure, Full Back
Seams	Sonic-Welded
Packaging	10/Bag, 10 Bags, 100/Case

This product is NOT made with natural rubber latex

*AATCC 42 - Water Resistance: Spray Impact Penetration Test
AATCC 127 - Water Resistance: Hydrostatic Pressure Test
According to AAMI Standard PB70:2012 Liquid Barrier Performance

For more information on these test methods, please visit aatcc.org and aami.org

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.

