

65-3112

Comfort-Plus™ Surgical Mask with Stretch Knit Ties


- Complies with ASTM F2100 Level 1 requirements
- Soft sewn stretch knit ties provide secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	Comfort-Plus™ Surgical Mask with Stretch Knit Ties
Color	Blue
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.6
Filtration - PFE	≥99.8
Delta P	≤2.6
Flammability	Class 1
Attachment Design	Vertical Stretch Knit Ties
Mask Construction	Sewn
Packaging	50/Box, 5/Boxes, 250/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 80
BFE (%) @ 3.0 µm: ≥ 95%
PFE (%) @ 0.1 µm: ≥ 95%
ΔP (mm H₂O/cm²): < 4
Flame Spread: Class 1