

65-3122

Comfort-Plus[™] Surgical Mask with Stretch Knit Ties


- Complies with ASTM F2100 Level 1 Requirements
- Soft sewn stretch knit ties provide secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	Comfort-Plus [™] Surgical Mask w/Stretch Knit Ties
Color	Green Diamond
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.6
Filtration - PFE	≥99.8
Delta P	≤2.6
Flammability	Class 1
Attachment Design	Vertical Stretch Knit Ties
Mask Construction	Sewn
Packaging	50/Bx, 5/Boxes, 250/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 μm:	≥ 95%
PFE (%) @ 0.1 μm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1