

65-3130

Sensitive Skin Plus Surgical Mask


- Complies with ASTM F2100 Level 1 Requirements
- Color-free face mask with specially designed inner facing improves comfort
- Soft sewn ties and bindings promote secure comfortable fit
- Integrated layers work together to maximize breathability



Specifications

Product	Sensitive Skin Plus Surgical Mask
Color	White
Fluid Resistance	80 mm Hg
Filtration - BFE	≥98.3
Filtration - PFE	≥98.1
Delta P	≤2.7
Flammability	Class 1
Attachment Design	Vertical Soft Poly Ties
Mask Construction	Sewn
Packaging	50/Box, 6/Boxes, 300/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 μm:	≥ 95%
PFE (%) @ 0.1 μm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1