

65-3312

DualGard™ High Filtration Surgical Mask


- Complies with ASTM F2100 Level 3 Requirements
- Soft sewn stretch knit ties provide secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	DualGard™ High Filtration Surgical Mask
Color	Blue Diamond
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99.8
Filtration - PFE	≥99.0
Delta P	≤3.2
Flammability	Class 1
Attachment Design	Vertical Stretch Knit Ties
Mask Construction	Sewn
Packaging	50/Box, 5/Boxes, 250/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 160
BFE (%) @ 3.0 μm: ≥ 98%
PFE (%) @ 0.1 μm: ≥ 98%
ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1