

65-3321

Tape Fog-Shield® Surgical Mask


- Complies with ASTM F2100 Level 1 Requirements
- Fog-Shield® anti-fog tape inhibits eyewear fogging improving vision
- Soft sewn ties and bindings promote secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	Tape Fog-Shield® Surgical Mask
Color	Green Diamond
Fluid Resistance	80 mm Hg
Filtration - BFE	≥99.6
Filtration - PFE	≥99.8
Delta P	≤2.6
Flammability	Class 1
Attachment Design	Vertical Soft Poly Ties
Mask Construction	Sewn
Packaging	50/Bx, 5/Boxes, 250/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 µm:	≥ 95%
PFE (%) @ 0.1 µm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1