

65-3344

DualGard™ Surgical Mask with Extended Shield


- Complies with ASTM F2100 Level 3 Requirements
- Snap-on coated film shield with proprietary anti-fog treatment provides enhanced clarity
- Snap-on extended shield can be removed from mask if needed
- Soft sewn ties and bindings promote secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	DualGard™ Surgical Mask with Extended Shield
Color	Blue Diamond
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99.9
Filtration - PFE	≥99.2
Delta P	≤2.8
Flammability	Class 1
Attachment Design	Vertical Stretch Knit Ties
Mask Construction	Sewn
Packaging	25/Box, 4/Boxes, 100/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg): 160
BFE (%) @ 3.0 μm: ≥ 98%
PFE (%) @ 0.1 μm: ≥ 98%
ΔP (mm H₂O/cm²): < 5
Flame Spread: Class 1