

65-3348

DualGard™ Foam Shield™ Surgical Mask with Stretch Knit Ties


- Complies with ASTM F2100 Level 3 Requirements
- Foam Shield™ provides added comfort with a foam strip that keeps shield off forehead
- Snap-on anti-fog Foam Shield can be removed from mask if needed
- Soft sewn stretch knit ties provide secure comfortable fit
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability



Specifications

Product	DualGard™ Foam Shield™ Surgical Mask with Stretch Knit Ties
Color	Blue Diamond
Fluid Resistance	160 mm Hg
Filtration - BFE	≥99.9
Filtration - PFE	≥99.2
Delta P	≤2.8
Flammability	Class 1
Attachment Design	Vertical Stretch Knit Ties
Mask Construction	Sewn
Packaging	25/Box, 4/Boxes, 100/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	160
BFE (%) @ 3.0 µm:	≥ 98%
PFE (%) @ 0.1 µm:	≥ 98%
ΔP (mm H ₂ O/cm ²):	< 5
Flame Spread:	Class 1