

15206

Lite & Cool Surgical Mask with Anti-Glare Shield


- Complies with ASTM F2100 Level 1 Requirements
- Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity
- Specially designed inner facing improves comfort
- Integrated layers work together to maximize breathability
- Soft longer ties promote secure comfortable fit



Specifications

Product	Lite & Cool Surgical Mask with Anti-Glare Shield
Color	Blue
Fluid Resistance	80 mm Hg
Filtration - BFE	≥98.3
Filtration - PFE	≥98.1
Delta P	≤2.7
Flammability	Class 1
Attachment Design	Vertical Ties
Mask Construction	Sonic Welded
Packaging	25/Box, 4/Boxes, 100/Cs

This product is NOT made with natural rubber latex



Complies with ASTM F2100

Fluid Resistance (mmHg):	80
BFE (%) @ 3.0 µm:	≥ 95%
PFE (%) @ 0.1 µm:	≥ 95%
ΔP (mm H ₂ O/cm ²):	< 4
Flame Spread:	Class 1