

51116

Isolation Gown

Full Coverage

- Polycoated material passes AATCC 42 and AATCC 127* testing.
- Material is low linting
- Knit cuffs are comfortable and facilitate glove donning.
- Full coverage with overlapping back panels.



Specifications

Product	Isolation Gown, Full Coverage
Size	X-Large
Dimensions	Length 47", Width 55"
Color	Yellow
Wrist	Knit Cuffs
Material	Polycoated
Design	Secure Neck Ties, Full Back
Seams	Sonic-Welded
Packaging	10/Bag, 5/Bags, 50/Case

This product is NOT made with natural rubber latex

***AATCC 42 - Spray Impact \leq 1.0g**
AATCC 127 - Hydrostatic Pressure \geq 50cm

For more information on these test methods, please visit aatcc.org and aami.org

Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards.