

Fluid Resistant Surgical Masks

ASTM F2100, Level 1 (80 mmHg)

Precept ASTM F2100 Level 1 face mask are designed for lower risk fluid environments, providing good protection with filtration efficiencies (BFE and PFE) exceeding ≥97%, and all masks meet Class 1 flame resistance requirements. Precept Level 1 face mask are available in a variety of styles including surgical and procedure with and without extended shields. All masks are sonic welded and have an enclosed malleable nosepieces to ensure a snug, yet comfortable fit.



15200



15201











65-3010



65-3110

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
--------	-------------	----------	-------	-----------

Fluid Resistant Surgical Face Masks

15200 	Lite & Cool Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Specially designed inner facing improves comfort Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	Blue 	50/Box, 6/Boxes, 300/Cs
15201 	Standard Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Specially designed inner facing improves comfort Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	Blue 	50/Box, 6/Boxes, 300/Cs
65-3010 	Comfort-Cool™ Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Soft sewn ties and bindings promote secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Blue 	50/Box, 6/Boxes, 300/Cs
65-3110 	Comfort-Plus™ Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Soft sewn ties and bindings promote secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Blue 	50/Box, 6/Boxes, 300/Cs