

Complies with ASTM F2100

Fluid Resistance (mmHg): 80

BFE (%) @ 3.0 µm: ≥ 95%

PFE (%) @ 0.1 µm: ≥ 95%

ΔP (mm H₂O/cm²): < 4

Flame Spread: Class 1

Fluid Resistant Specialty Surgical Masks

ASTM F2100, Level 1 (80 mmHg)

Precept ASTM F2100 Level 1 face mask are designed for lower risk fluid environments, providing good protection with filtration efficiencies (BFE and PFE) exceeding ≥97%, and all masks meet Class 1 flame resistance requirements. Precept Level 1 face mask are available in a variety of styles including surgical and procedure with and without extended shields. All masks are sonic welded and have an enclosed malleable nosepieces to ensure a snug, yet comfortable fit.



15210



15211




15212

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
--------	-------------	----------	-------	-----------

Fluid Resistant Specialty Surgical Face Masks

15210 	Foam Anti-Fog Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Fog-Shield® anti-fog foam inhibits eyewear fogging improving vision Specially designed inner facing improves comfort Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	Blue 	50/Box, 6/Boxes, 300/Cs
15211 	Film Anti-Fog Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Anti-fog film inhibits eyewear fogging improving vision Specially designed inner facing improves comfort Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	Green 	50/Box, 6/Boxes, 300/Cs
15212 	Tape Anti-Fog Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Anti-fog tape inhibits eyewear fogging improving vision Specially designed inner facing improves comfort Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	Green 	50/Box, 6/Boxes, 300/Cs



Complies with ASTM F2100

Fluid Resistance (mmHg): 80

BFE (%) @ 3.0 µm: ≥ 95%

PFE (%) @ 0.1 µm: ≥ 95%

ΔP (mm H₂O/cm²): < 4

Flame Spread: Class 1

Fluid Resistant Specialty Surgical Masks

ASTM F2100, Level 1 (80 mmHg)



15215



65-3112



65-3120

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
Fluid Resistant Specialty Surgical Face Masks				
15215 	Sensitive Skin Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Color-free face mask with specially designed inner facing improves comfort Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	White 	50/Box, 6/Boxes, 300/Cs
15216 	Sensitive Skin Foam Anti-Fog Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Color-free face mask with specially designed inner facing improves comfort Anti-fog foam inhibits eyewear fogging improving vision Integrated layers work together to maximize breathability Soft longer ties promote secure comfortable fit 	White 	50/Box, 6/Boxes, 300/Cs
65-3112 	Comfort-Plus™ Surgical Mask with Stretch Knit Ties	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Soft sewn stretch knit ties provide secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Blue 	50/Box, 5/Boxes, 250/Cs
65-3120 	Comfort-Plus™ Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Soft sewn ties and bindings promote secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Green 	50/Box, 6/Boxes, 300/Cs

Fluid Resistant Specialty Surgical Masks

ASTM F2100, Level 1 (80 mmHg)









65-3122




65-3130



65-3320

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
Fluid Resistant Specialty Surgical Face Masks				
65-3122 	Comfort-Plus™ Surgical Mask with Stretch Knit Ties	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Soft sewn stretch knit ties provide secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Green 	50/Box, 5/Boxes, 250/Cs
65-3130 	Sensitive Skin Plus Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Color-free face mask with specially designed inner facing improves comfort Soft sewn ties and bindings promote secure comfortable fit Integrated layers work together to maximize breathability 	White 	50/Box, 6/Boxes, 300/Cs
65-3320 	Foam Fog-Shield® Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Fog-Shield® anti-fog foam inhibits eyewear fogging improving vision Strong polyester ties provide secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Green 	25/Box, 6/Boxes, 150/Cs



Complies with ASTM F2100

Fluid Resistance (mmHg): 80

BFE (%) @ 3.0 µm: ≥ 95%

PFE (%) @ 0.1 µm: ≥ 95%

ΔP (mm H₂O/cm²): < 4

Flame Spread: Class 1

Fluid Resistant Specialty Surgical Masks





ASTM F2100, Level 1 (80 mmHg)



65-3321



65-3322

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
65-3321 	Tape Fog-Shield® Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Fog-Shield® anti-fog tape inhibits eyewear fogging improving vision Soft sewn ties and bindings promote secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Green 	50/Box, 5/Boxes, 250/Cs
65-3322 	Tape Fog-Shield® Surgical Mask	<ul style="list-style-type: none"> Complies with ASTM F2100 Level 1 Requirements Fog-Shield® anti-fog tape inhibits eyewear fogging improving vision Strong polyester ties provide secure comfortable fit Specially designed inner facing improves comfort Integrated layers work together to maximize breathability 	Green 	50/Box, 5/Boxes, 250/Cs