



**Complies with ASTM F2100**

Fluid Resistance (mmHg): 120

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1

## Fluid Resistant Surgical Masks

ASTM F2100, Level 2 (120 mmHg)

**Precept** ASTM F2100 Level 2 face mask are designed for moderate to high risk fluid environments, providing excellent protection with filtration efficiencies (BFE and PFE) exceeding ≥98%, and all masks meet Class 1 flame resistance requirements. Precept Level 2 face mask are available in a variety of styles including surgical and procedure with and without extended shields. All masks are sonic welded and have an enclosed malleable nosepieces to ensure a snug, yet comfortable fit.



15251









15820



15920

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
--------	-------------	----------	-------	-----------

### Fluid Resistant Surgical Face Masks

15251 	FluidGuard® 120 Horizontal Tie Surgical Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 2 Requirements</li> <li>Wider design provides more facial coverage than standard sized masks</li> <li>Soft longer horizontal ties promote secure comfortable fit</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> </ul>	Blue with Pattern 	50/Box, 12/Boxes, 600/Cs
15820 	FluidGuard® 120 FluidPouch™ Surgical Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 2 Requirements</li> <li>Pouch style mask provides off the face design giving more room to breathe</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft longer ties promote secure comfortable fit</li> </ul>	Blue with Pattern 	50/Box, 6/Boxes, 300/Cs
15920 	FluidGuard® 120 Anti-Fog Surgical Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 2 Requirements</li> <li>Anti-fog foam inhibits eyewear fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft longer ties promote secure comfortable fit</li> </ul>	Blue 	50/Box, 6/Boxes, 300/Cs



Complies with ASTM F2100

Fluid Resistance (mmHg): 120

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1

## Fluid Resistant Surgical Masks

ASTM F2100, Level 2 (120 mmHg)



15255



15830



15930

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
--------	-------------	----------	-------	-----------

### Fluid Resistant Surgical Face Masks with Shields

15255 	FluidGuard® 120 Horizontal Tie Surgical Mask with Anti-Glare Shield	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 2 Requirements</li> <li>Coated film shield with proprietary anti-fog, anti-glare treatment provides enhanced clarity</li> <li>Wider design provides more facial coverage than standard sized masks</li> <li>Soft longer horizontal ties promote secure comfortable fit</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> </ul>	Blue with Pattern 	25/Box, 4/Boxes, 100/Cs
15830 	FluidGuard® 120 FluidPouch™ Surgical Mask with Anti-Glare Shield	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 2 Requirements</li> <li>Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity</li> <li>Pouch style mask provides off the face design giving more room to breathe</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft longer ties promote secure comfortable fit</li> </ul>	Blue with Pattern 	25/Box, 4/Boxes, 100/Cs
15930 	FluidGuard® 120 Anti-Fog Surgical Mask with Extended Shield	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 2 Requirements</li> <li>Coated film shield with proprietary anti-fog treatment provides enhanced clarity</li> <li>Anti-fog foam inhibits fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft longer ties promote secure comfortable fit</li> </ul>	Blue 	25/Box, 4/Boxes, 100/Cs