

**Complies with ASTM F2100**

Fluid Resistance (mmHg): 160

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1

## Fluid Resistant Procedure Masks

ASTM F2100, Level 3 (160 mmHg)

**Precept** ASTM F2100 Level 3 face mask are designed for high risk fluid environments, providing excellent protection with filtration efficiencies (BFE and PFE) exceeding ≥98%, and all masks meet Class 1 flame resistance requirements. Precept Level 3 face mask are available in a variety of styles including surgical and procedure with and without extended shields. All masks are sonic welded and have an enclosed malleable nosepieces to ensure a snug, yet comfortable fit.



14400



15300



15301


ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
--------	-------------	----------	-------	-----------

### Fluid Resistant Procedure Face Masks

14400 	FluidGuard® 160 Procedure Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue 	50/Box, 10/Boxes, 500/Cs
14401 	FluidGuard® 160 Anti-Fog Procedure Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Anti-fog foam inhibits eyewear fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue 	50/Box, 10/Boxes, 500/Cs
15300 	FluidGuard® 160 Procedure Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue Diamond 	50/Box, 10/Boxes, 500/Cs
15301 	FluidGuard® 160 Anti-Fog Procedure Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Anti-fog foam inhibits eyewear fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue Diamond 	50/Box, 10/Boxes, 500/Cs

Precept products do NOT contain natural rubber latex

For more information or to order, contact Precept at **800.438.5827** or visit [www.preceptmed.com](http://www.preceptmed.com)



Complies with ASTM F2100

Fluid Resistance (mmHg): 160

BFE (%) @ 3.0 µm: ≥ 98%

PFE (%) @ 0.1 µm: ≥ 98%

ΔP (mm H<sub>2</sub>O/cm<sup>2</sup>): < 5

Flame Spread: Class 1

## Fluid Resistant Procedure Masks

ASTM F2100, Level 3 (160 mmHg)



14410



15311



65-3347

ITEM #	DESCRIPTION	FEATURES	COLOR	PACKAGING
--------	-------------	----------	-------	-----------

### Fluid Resistant Procedure Face Masks with Shields

14410 	FluidGuard® 160 Anti-Fog Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Coated film shield with proprietary anti-fog treatment provides enhanced clarity</li> <li>Anti-fog foam inhibits fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue 	25/Box, 4/Boxes, 100/Cs
15310 	FluidGuard® 160 Anti-Fog Procedure Mask with Extended Shield	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Coated film shield with proprietary anti-fog treatment provides enhanced clarity</li> <li>Anti-fog foam inhibits fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue Diamond 	25/Box, 4/Boxes, 100/Cs
15311 	FluidGuard® 160 Anti-Fog Procedure Mask with Anti-Glare Shield	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Coated film shield with proprietary anti-glare, anti-fog treatment provides enhanced clarity</li> <li>Anti-fog foam inhibits fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue Diamond 	25/Box, 4/Boxes, 100/Cs
65-3347 	DualGuard™ Foam Shield™ Anti-Fog Procedure Mask	<ul style="list-style-type: none"> <li>Complies with ASTM F2100 Level 3 Requirements</li> <li>Foam Shield™ provides added comfort with a foam strip that keeps shield off forehead</li> <li>Snap-on anti-fog Foam Shield can be removed from mask if needed</li> <li>Anti-fog foam inhibits fogging improving vision</li> <li>Specially designed inner facing improves comfort</li> <li>Integrated layers work together to maximize breathability</li> <li>Soft comfortable earloops provide easier donning</li> </ul>	Blue 	25/Box, 4/Boxes, 100/Cs

Precept products do NOT contain natural rubber latex

For more information or to order, contact Precept at **800.438.5827** or visit [www.preceptmed.com](http://www.preceptmed.com)