



51196 Isolation Gown – full coverage AAMI level 2

Intended Use

Intended to be used in environments where low to moderate fluid exposure is a risk.

Key Features

- Multi-layer SMS material[®] delivers AAMI Level 2 protection and comfort
- Full coverage with overlapping panels
- Knit cuffs are comfortable and facilitate glove donning
- Secure neck ties and sonic-welded seams
- Universal, one size fits most
- Made in North America
- Not made with natural rubber latex



Barrier Protection AAMI PB70:2012

Test Method	Test Result	Anticipated Risk of Exposure
AATCC 42 Spray Impact	≤ 1.0g	Low
AATCC 127 Hydrostatic Pressure	≥ 20cm	Low

Isolation Gown – full coverage AAMI level 2

Intended to be used in environments where low to moderate fluid exposure is a risk.

Technical Specifications

Part No.	51196
Description	Isolation Gown, Full Coverage
Color	Blue
Size	Universal
Dimensions	Length 48", Width 60"
Wrist	Knit Cuffs
Material	Multi-Layer SMS
Design	Secure Neck Ties, Full Coverage
Seams	Sonic-Welded
Qty	10/Bag, 10/Bags, 100/Case



ⓐAATCC 42 – Water Resistance: Spray Impact Penetration Test / AATCC 127 – Water Resistance: Hydrostatic Pressure Test For more information on these test methods, please visit aatcc.org and aami.org Precept makes no representation or warranty that its products will provide absolute protection from contact or exposure to contaminants or blood, or eliminate the risk of contracting diseases of any nature. OSHA requires the employer to evaluate the task and type of exposure anticipated and select the appropriate personal protective equipment (PPE) that reduces the employee's exposure to potential hazards. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2020 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.
PA0331 rev 1 2001 25-SEPT-2020 ENG – US

**Precept Medical Products, Inc. an
Aspen Surgical Products, Inc. Company**
phone +1 800.438.5827
preceptmed.com

